(43) International Publication Date 13 January 2005 (13.01.2005)

PCT

(10) International Publication Number WO 2005/002472 A1

- (51) International Patent Classification⁷: A611, 27/58
- A61F 2/04.
- (21) International Application Number:

PCT/SE2004/000879

- (22) International Filing Date:
- 7 June 2004 (07.06.2004)
- (25) Flüng Language:

English

(26) Publication Language:

English

(30) Priority Data: 0301974-2


3 July 2003 (03.07.2003) SE

- (71) Applicant (for all designated States except US): ASTRA
 TECH AB [SE/SE]; Aminogatan 1, S-431 53 Mölndal
 (SE).
- (72) laventors; and
- (75) Inventors/Applicants (for t/S only): OLAUSSON, Susanne [SE/SE]; Holmenbacken 19B, S-443 38 Lerum (SE). WIBERG, Mikael [SE/SE]: Prästgården, S-870 33 Docksta (SE). RINDLAV WESTLING, Åss [SE/SE]; Vallis Parkväg 42, S-437 32 Lindorne (SE).

- (74) Agent: AWAPATENT AB; P.O. Box 11394, S-404 28 Stockholm (SE).
- (81) Designated States (unless otherwise indicated, for every kind of national protection available): AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BW, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, BG, BS, Fl, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JF, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NA, NI, NO, NZ, OM, PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, SY, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW.
- (84) Designated States (unless otherwise indicated, for every kind of regional protection available): ARIPO (BW, GH, GM, KE, LS, MW, MZ, NA, SD, SL, SZ, TZ, UG, ZM, ZWI, Eurasian (AM, AZ, BY, KG, KZ; MD, RU, TJ, TM, European (AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES, FR, GB, GR, HU, IE, IT, LU, MC, NL, PL, PT, RO, SE, SI, SK, TR, OAPI (BF, BJ, CF, CG, Cl, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG).

[Continued on next page]

(54) Title: A DEVICE FOR PROMOTING REGENERATION OF AN INJURED NERVE, A KIT AND A BIODEGRADABLE SHEET FOR PREPARING SUCH A DEVICE



(57) Abstract: The invention relates to promotion of a process for regeneration of an injured nerve using a plurality of guiding means, preferably a plurality of guiding fibres, presenting an in vivo biodegradability being such that at least a majority of said guiding means becomes essentially disintegrated by degradation (and/or dissolution) during a pre-contact period extending from application of the device at the injured nerve up to a first occurrence of a re-established (regenerated) contact between the ends of the injured nerve. During a post-contact period extending from the end of the pre-contact period and up to the end of the regeneration process, the disintegrated guiding means will provide no substantial axon growth guiding function and no substantial axon growth blocking effect.

BEST AVAILABLE COPY